

TESLAFormer™

4 CHANNELS

**BUILD MUSCLE  
BURN FAT**  
*with*  
FUNCTIONAL  
MAGNETIC  
STIMULATION



# TESLAFormer™

## HOW IT WORKS

TESLAFormer™ is a more concentrated, targeted and effective device for body forming. The magnetic field propagates through clothes and skin inside your body. It is used for muscle strengthening on areas of glutes, abdomen, hamstring area, arms and pelvic floor using the new technology of FMS (Functional Magnetic Stimulation).

TESLAFormer™ will provide a pulsed magnetic field that will contract the muscles through its effect on motor nerves of the muscles on the target site. It has been scientifically found that the technology of FMS will result in increased muscle mass after multiple treatments on the same target location.

## ADVANTAGES

- ✓ Growing and strengthening muscle
- ✓ Non-invasive treatment
- ✓ No downtime
- ✓ More than 50,000 contractions in a 30 minute session
- ✓ Electro magnetic energy
- ✓ No operator needed
- ✓ No preparation needed

## KEY BENEFITS

- ✓ To form your body
- ✓ Grow muscle mass
- ✓ Easier fat burning
- ✓ Stay in shape
- ✓ Butt lifting
- ✓ Assist posture



# TESLAChair™

## FUNCTIONAL MAGNETIC STIMULATION

*A breakthrough treatment for incontinence and confidence. This unique technology revolutionises the women's intimate health and wellness category by providing those suffering from incontinence with a **completely non-invasive option**.*



## CHAIR FOR PELVIC AREA FOR THE REPAIR OF INCONTINENCE

The TESLAChair™ was developed for treating weakened pelvic floor muscles and urinary incontinence using the special Functional Magnetic Stimulation (FMS) technology. The pulsed magnetic field generated by the device causes the muscles of the pelvic floor to contract. This therapy is suitable for all kinds of urinary and fecal incontinence. FMS also helps patients with sexual dysfunction and chronic pain in the pelvic region and lumbar spine. The magnetic stimulation results in increased strength and endurance of the sphincter and pelvic floor muscles. The therapy also has a positive effect on sexual performance, erections, the ability to achieve orgasm and the active contracting of groin and gluteal muscles.


## CHAIR FOR LOWER BACK

The new chair for treatment of the lower back area is a painless, non-invasive and effective therapy for problems in the lower back area. The treatment stimulates muscles, and accelerates circulation without direct contact with the skin. The therapy procedure is simple, and the patient does not require special preparation.


During the therapy the patient simply sits in the chair. Stimulation of the lower back can be accurately targeted by changing the position of the applicator in the chair backrest.

## TESLAChair™ INCONTINENCE & PELVIC FLOOR


- Urge urinary incontinence
- Stress urinary incontinence
- Mixed urinary incontinence
- Incontinence after radical prostatectomy
- Fecal incontinence
- Regeneration after childbirth
- Milder forms of drooping uterus and bladder (cistokela, descensus uterusa)
- Rehabilitation after pelvic floor surgery
- Chronic prostate inflammation
- Erectile dysfunction



Chair for pelvic area (1600632)



Chair for lumbar area (1600638)



Chair for pelvic & lumbar area (1600635)

**STOP THE URGE  
BUILD PELVIC STRENGTH**

# TESLAFormer™

## TECHNICAL DATA

MAGNETIC FIELD STRENGTH	Up to 3 Tesla
THERAPY FREQUENCY	1 to 160 Hz
NUMBER OF CHANNELS	4
INTERFACE	Colour touchscreen
MODULATION	Amplitude or frequency modulation
ACTIVE TIME	1 s - 20 s
PAUSE TIME	0 s - 240 s
MAX NO. OF STEPS PER PROGRAM	8
MAX DURATION OF EACH STEP	1 s - 20 min
THERAPY DURATION	1 min - 60 min
MANUAL PROGRAMS	Yes
PRE-SET PROGRAMS	Yes
PRE-SET PROTOCOLS	Yes
BODY PART PROGRAMS	Yes
VOLTAGE, POWER	AC240V 50-60 Hz, 1500 VA
SIZE: TESLACare / TESLACare prestige	40x51x26 / 47x60x27 cm (WxLxH)
WEIGHT: TESLACare / TESLACare prestige	9.8 / 10.7 kg
MDD 93/42/EEC classification	Class IIa/ARTG Inclusion No. 317701



Chair for Pelvic & Back 2 Prestige (1600649)