

3-in-1  
super-facial

## GENEO BENEFITS

### Immediate results that stay

Effective formulations and treatment kits that provide immediate and long-lasting results, for revitalizing, purifying and rejuvenating the skin, and achieving overall improvement of skin texture and appearance.

### Personalized and versatile

A single platform containing 3 powerful technologies that can be used in synergy to answer specific needs and desires.

### Safe and pleasant

Tested and validated technologies, designed to stimulate the body's natural healing mechanisms, with minimal discomfort and downtime.

geneo

### About Pollogen™

Pollogen Ltd. a company of Lumenis, is a global leader in the medical and home-use aesthetics devices markets, providing innovative, safe, and effective solutions that enhance the practice of medical aesthetics professionals, as well as serving private aesthetic consumers. The company offers products that support a range of skin treatments, based on its proprietary technologies which include TriPollar, Voulderm, Trifractional, OxyGeneo and DMA. Pollogen provides solutions for professionals and consumers in more than 60 countries, addressing the specific needs and the ever-changing standards of professionals and consumers in the dynamic beauty industry.

#### GENEO TECHNICAL SPECIFICATIONS

**WEIGHT** 4.4 Kg

**OVERALL DIMENSIONS** 380MM x 382MM x 282MM (H x W x D)

**INPUT VOLTAGE OF POWER SUPPLY** 100-240V~50/60Hz

**INPUT POWER RATING OF POWER SUPPLY** 12VDC, 2.0A

#### TRIPOLLAR APPLICATOR

**MAXIMUM POWER OUTPUT** Up to 15 Watts @ 200 Ohm

**OUTPUT FREQUENCY** 1MHz

**OUTPUT POWER CONTROL** Pulse width modulation (PWM)

### About Lumenis®

Lumenis is a global leader in the field of minimally-invasive clinical solutions for the Surgical, Ophthalmology and Aesthetic markets, and a world renowned expert in developing and commercializing innovative energy-based technologies, including Laser, Intense Pulsed Light (IPL) and Radio-Frequency (RF).

#### OXYGENEO APPLICATOR

**POWER OUTPUT** Vibration - Frequency 5 – 30Hz

#### ULTRASOUND APPLICATOR

**MAXIMUM POWER OUTPUT** Up to 1 Watts /cm2

**OUTPUT FREQUENCY** 1MHz

**OUTPUT POWER CONTROL** Pulse width modulation (PWM)

© 2019 Pollogen Ltd. All rights reserved.

Trademarks and logos of geneo, TriPollar, OxyGeneo and Pollogen are all trademarks of Pollogen Ltd.

11600901 | REV1 | APRIL 2019



geneo



breathe in,  
shine out

## INTRODUCING GENE0

Geneo is the ultimate facial treatment platform for achieving a younger, fresher and naturally better-looking skin. It offers a unique, three-in-one facial treatment that works both on the outside and the inside of the skin, combining clinically proven technologies: Tripollar RF, Oxygeneo and Ultrasound. With these technologies we encourage the natural processes of the body to renew and replenish the skin.



## HOW DOES IT WORK?

Geneo consists of 3 applicators offering the singular benefits of 3 potent technologies – OxyGeneo, TriPollar RF and Ultrasound. Each technology can be used separately to enjoy its benefits, or combined for a full synergetic effect.

### Skin Replenishment

Collagen and Elastin are major building blocks of the skin, and their depletion results in fine lines and wrinkles. With Geneo Tripollar RF applicator users enjoy the opportunity to encourage the natural remodeling and replenishment processes of Collagen and Elastin fibers. This treatment applies gentle heat to the layers of the skin, resulting in Improving the appearance of wrinkles and fine lines, and visible skin firmness and lifting effect.


### Skin Revitalization

The revolutionary Geneo treatment offers the unique triple effect of natural oxygenation, exfoliation and absorption of active ingredients. By applying a gentle massage with the capsule and the gel, a reaction which utilizes natural physiological processes is created, naturally elevating oxygen levels in the skin. The treatment is finished with the dedicated serum, introducing beneficial ingredients to the skin in its highly receptive state. This pleasant, refreshing treatment results in immediate revitalization of the skin from the very first treatment, which is further

strengthened and maintained throughout the recommended series of treatments.

### Skin Enhancement

By introducing Ultrasound waves to the skin, Geneo enhances the effect of replenishment and revitalization. The ultrasound waves improve blood flow and assist with better absorption of active ingredients. This safe and painless treatment results in visible improvement of skin texture and tone.



## GENEO TECHNOLOGIES

### OXYGENEO®



Refreshing and improving skin texture and tone with the 3-in-1 concept: by exploiting the physiological processes responsible for increasing in oxygen levels in the skin, gently exfoliating the skin and preparing it for the absorption of active ingredients.

### TRIPOLLAR® RF



Triggering the renewal of the skin by stimulating production of Collagen and Elastin fibers, via the delivery of focused RF energy in a safe and effective manner, inducing controlled heating of the dermis.

### ULTRASOUND



Improving blood flow and absorption of active ingredients by applying ultrasound waves to the skin, creating vibration and re-arrangement of the epidermal cells in order to increase the permeability of the stratum corneum and improve skin nutrition.